



groL Antibody

Product Code	CSB-PA330104HA01DRZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31681
Immunogen	Recombinant Chlamydia pneumoniae 60 kDa chaperonin protein (1-214AA)
Raised In	Rabbit
Species Reactivity	Chlamydia pneumoniae
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60 kDa chaperonin (GroEL protein) (Protein Cpn60), groL, groEL mopA
Immunogen Species	Chlamydia pneumoniae (Chlamydophila pneumoniae)
Research Area	Others
Target Names	groL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.