




pbpE Antibody

Product Code	CSB-PA333755HA01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32959
Immunogen	Recombinant Bacillus subtilis Penicillin-binding protein 4* protein (1-451AA)
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Penicillin-binding protein 4* (PBP 4*) (PBP 4A) (Penicillin-binding protein E), pbpE
Immunogen Species	Bacillus subtilis
Research Area	Others
Target Names	pbpE

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: pbpE antibody at 3μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 52 kDa

Observed band size: 52 kDa