



pbpE Antibody, HRP conjugated

Product Code	CSB-PA333755HB01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32959
Immunogen	Recombinant Bacillus subtilis Penicillin-binding protein 4* protein (1-451AA)
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Penicillin-binding protein 4* (PBP 4*) (PBP 4A) (Penicillin-binding protein E), pbpE
Immunogen Species	Bacillus subtilis (strain 168)
Research Area	Others
Target Names	pbpE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.