



PPX1 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA334459HC01SVG |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P38698 |
| Immunogen | Recombinant Saccharomyces cerevisiae Exopolyphosphatase protein (1-397AA) |
| Raised In | Rabbit |
| Species Reactivity | Saccharomyces cerevisiae |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Exopolyphosphatase (ExopolyPase) (EC 3.6.1.11) (Metaphosphatase), PPX1 |
| Immunogen Species | Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast) |
| Research Area | Others |
| Target Names | PPX1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |