



BETVII Antibody

Product Code	CSB-PA335151HA01BSS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25816
Immunogen	Recombinant Betula pendula Profilin-1 protein (2-133AA)
Raised In	Rabbit
Species Reactivity	Betula pendula
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Profilin-1 (Allergen Bet v II) (Pollen allergen Bet v 2) (allergen Bet v 2), BETVII
Immunogen Species	Betula pendula (European white birch) (Betula verrucosa)
Research Area	Others
Target Names	BETVII
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.