




PHLP II Antibody

Product Code	CSB-PA337346HA01EUQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43214
Immunogen	Recombinant Phleum pratense Pollen allergen Phl p 2 protein (27-122AA)
Raised In	Rabbit
Species Reactivity	Phleum pratense
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pollen allergen Phl p 2 (Allergen Phl p II) (allergen Phl p 2), PHLPII
Immunogen Species	Phleum pratense
Research Area	Others
Target Names	PHLP II

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: PHLPII antibody at 2.8µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 27 kDa

Observed band size: 27 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.