



pepF Antibody, HRP conjugated

Product Code	CSB-PA346917LB01MLW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54125
Immunogen	Recombinant Mycoplasma pneumoniae Oligoendopeptidase F homolog protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Mycoplasma pneumoniae
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oligoendopeptidase F homolog (EC 3.4.24.-), pepF
Immunogen Species	Mycoplasma pneumoniae (strain ATCC 29342 / M129)
Research Area	Others
Target Names	pepF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.