



ompT Antibody, Biotin conjugated

Product Code	CSB-PA348292LD01EOD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58603
Immunogen	Recombinant Escherichia coli O157:H7 Protease 7 protein (21-317AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli O157:H7
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protease 7 (EC 3.4.23.49) (OmpT) (Outer membrane protein 3B) (Protease A) (Protease VII), ompT
Immunogen Species	Escherichia coli O157:H7
Research Area	Others
Target Names	ompT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.