



ompK Antibody, HRP conjugated

Product Code	CSB-PA349512LB01VFE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P59570
Immunogen	Recombinant Vibrio parahaemolyticus serotype O3:K6 Outer membrane protein OmpK protein (21-266AA)
Raised In	Rabbit
Species Reactivity	Vibrio parahaemolyticus serotype O3:K6
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane protein OmpK, ompK
Immunogen Species	Vibrio parahaemolyticus serotype O3:K6 (strain RIMD 2210633)
Research Area	Others
Target Names	ompK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.