




ald1 Antibody

Product Code	CSB-PA350554HA01SKX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63478
Immunogen	Recombinant Staphylococcus aureus Alanine dehydrogenase 1 protein (1-372AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alanine dehydrogenase 1 (EC 1.4.1.1), ald1
Immunogen Species	Staphylococcus aureus
Research Area	Others
Target Names	ald1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: ald1 antibody at 2.4μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 57 kDa

observed band size: 57 kDa