




eno Antibody

Product Code	CSB-PA351974LA01SKX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P64078
Immunogen	Recombinant Staphylococcus aureus Enolase protein (1-434AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Enolase (EC 4.2.1.11) (2-phospho-D-glycerate hydro-lyase) (2-phosphoglycerate dehydratase), eno
Immunogen Species	Staphylococcus aureus
Research Area	Others
Target Names	eno

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: eno antibody at 3.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 64 kDa

Observed band size: 64 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.