



# eno Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA351974LC01SKX
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P64078
<b>Immunogen</b>	Recombinant Staphylococcus aureus Enolase protein (1-434AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Staphylococcus aureus
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Enolase (EC 4.2.1.11) (2-phospho-D-glycerate hydro-lyase) (2-phosphoglycerate dehydratase), eno
<b>Immunogen Species</b>	Staphylococcus aureus (strain Mu50 / ATCC 700699)
<b>Research Area</b>	Others
<b>Target Names</b>	eno
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.