




rpIL Antibody

Product Code	CSB-PA354049HA01SKX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P66061
Immunogen	Recombinant Staphylococcus aureus 50S ribosomal protein L7/L12 protein (1-122AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	50S ribosomal protein L7/L12, rpIL
Immunogen Species	Staphylococcus aureus
Research Area	Others
Target Names	rpIL

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: rpIL antibody at 3.25µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.