



rpIL Antibody, Biotin conjugated

Product Code	CSB-PA354049HD01SKX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P66061
Immunogen	Recombinant Staphylococcus aureus 50S ribosomal protein L7/L12 protein (1-122AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	50S ribosomal protein L7/L12, rpIL
Immunogen Species	Staphylococcus aureus (strain Mu50 / ATCC 700699)
Research Area	Others
Target Names	rpIL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.