



Beta-lytic metalloendopeptidase Antibody, HRP conjugated

Product Code	CSB-PA355563LB01LLB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00801
Immunogen	Recombinant Lysobacter enzymogenes Beta-lytic metalloendopeptidase protein (1-178AA)
Raised In	Rabbit
Species Reactivity	Lysobacter enzymogenes
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-lytic metalloendopeptidase (EC 3.4.24.32) (Beta-lytic protease)
Immunogen Species	Lysobacter enzymogenes
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.