



PKHD1 Antibody, Biotin conjugated

Product Code	CSB-PA357550LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08F94
Immunogen	Recombinant Human Fibrocystin protein (520-743AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibrocystin (Polycystic kidney and hepatic disease 1 protein) (Polyductin) (Tigmin), PKHD1, FCYT TIGM1
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	PKHD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.