




ompT Antibody

Product Code	CSB-PA357589LA01EVN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09169
Immunogen	Recombinant Escherichia coli Protease 7 protein (21-317AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:2000-1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protease 7 (EC 3.4.23.49) (OmpT) (Outer membrane protein 3B) (Protease A) (Protease VII), ompT
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	ompT

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: ompT antibody at 2.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 38 kDa

observed band size: 38 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.