



ompT Antibody, FITC conjugated

Product Code	CSB-PA357589LC01EVN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09169
Immunogen	Recombinant Escherichia coli Protease 7 protein (21-317AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protease 7 (EC 3.4.23.49) (OmpT) (Outer membrane protein 3B) (Protease A) (Protease VII), ompT
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	ompT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.