



hlgB Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA357862HA01SKY |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P0A075 |
| Immunogen | Recombinant Staphylococcus aureus Gamma-hemolysin component B protein (26-325AA) |
| Raised In | Rabbit |
| Species Reactivity | Staphylococcus aureus |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Gamma-hemolysin component B (H-gamma-1) (H-gamma-I), hlgB |
| Immunogen Species | Staphylococcus aureus (strain N315) |
| Research Area | Others |
| Target Names | hlgB |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |