



hlgB Antibody, Biotin conjugated

Product Code	CSB-PA357862HD01SKY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A075
Immunogen	Recombinant Staphylococcus aureus Gamma-hemolysin component B protein (26-325AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-hemolysin component B (H-gamma-1) (H-gamma-I), hlgB
Immunogen Species	Staphylococcus aureus (strain N315)
Research Area	Others
Target Names	hlgB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.