




Rubredoxin Antibody

Product Code	CSB-PA360387HA01CMA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00268
Immunogen	Recombinant Clostridium pasteurianum Rubredoxin protein (1-54AA)
Raised In	Rabbit
Species Reactivity	Clostridium pasteurianum
Tested Applications	ELISA, WB; Recommended dilution: WB:1:2000-1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rubredoxin (Rd)
Immunogen Species	Clostridium pasteurianum
Research Area	Others

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: Rubredoxin antibody at 2.8µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 22 kDa

observed band size: 22 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.