



Rubredoxin Antibody, Biotin conjugated

Product Code	CSB-PA360387HD01CMA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00268
Immunogen	Recombinant Clostridium pasteurianum Rubredoxin protein (1-54AA)
Raised In	Rabbit
Species Reactivity	Clostridium pasteurianum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rubredoxin (Rd)
Immunogen Species	Clostridium pasteurianum
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.