




RNU2 Antibody

Product Code	CSB-PA360432LA01UBD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00654
Immunogen	Recombinant Ustilago sphaerogena Ribonuclease U2 protein (1-114AA)
Raised In	Rabbit
Species Reactivity	Ustilago sphaerogena
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease U2 (RNase U2) (EC 3.1.27.4), RNU2
Immunogen Species	Ustilago sphaerogena
Research Area	Others
Target Names	RNU2

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: RNU2 antibody at 3.25µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 30 kDa

observed band size: 30 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.