



RNU2 Antibody, Biotin conjugated

Product Code	CSB-PA360432LD01UBD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00654
Immunogen	Recombinant Ustilago sphaerogena Ribonuclease U2 protein (1-114AA)
Raised In	Rabbit
Species Reactivity	Ustilago sphaerogena
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease U2 (RNase U2) (EC 3.1.27.4), RNU2
Immunogen Species	Ustilago sphaerogena (Smut fungus)
Research Area	Others
Target Names	RNU2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.