



gpd-1 Antibody, Biotin conjugated

Product Code	CSB-PA361431LD01CXY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04970
Immunogen	Recombinant Caenorhabditis elegans Glyceraldehyde-3-phosphate dehydrogenase 1 protein (1-341AA)
Raised In	Rabbit
Species Reactivity	Caenorhabditis elegans
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glyceraldehyde-3-phosphate dehydrogenase 1 (GAPDH-1) (EC 1.2.1.12), gpd-1
Immunogen Species	Caenorhabditis elegans
Research Area	Others
Target Names	gpd-1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.