



gB Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA362479LA01HJX |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P08666 |
| Immunogen | Recombinant Human herpesvirus 2 Envelope glycoprotein B protein (23-279AA) |
| Raised In | Rabbit |
| Species Reactivity | Human herpesvirus 2 |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Envelope glycoprotein B (gB), gB |
| Immunogen Species | Human herpesvirus 2 (strain HG52) (HHV-2) (Human herpes simplex virus 2) |
| Research Area | Others |
| Target Names | gB |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |