



gB Antibody, HRP conjugated

Product Code	CSB-PA362479LB01HJX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08666
Immunogen	Recombinant Human herpesvirus 2 Envelope glycoprotein B protein (23-279AA)
Raised In	Rabbit
Species Reactivity	Human herpesvirus 2
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Envelope glycoprotein B (gB), gB
Immunogen Species	Human herpesvirus 2 (strain HG52) (HHV-2) (Human herpes simplex virus 2)
Research Area	Others
Target Names	gB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.