




# ptsH Antibody

<b>Product Code</b>	CSB-PA362538HA01BRJ
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P08877
<b>Immunogen</b>	Recombinant Bacillus subtilis Phosphocarrier protein HPr protein (2-88AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Bacillus subtilis
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:2000-1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Phosphocarrier protein HPr (Histidine-containing protein), ptsH
<b>Immunogen Species</b>	Bacillus subtilis
<b>Research Area</b>	Others
<b>Target Names</b>	ptsH

## Image



### Western Blot

Positive WB detected in Recombinant protein

All lanes: ptsH antibody at 2.8µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 26 kDa

Observed band size: 26 kDa