




Arginine esterase Antibody

Product Code	CSB-PA362734LA01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09582
Immunogen	Recombinant Dog Arginine esterase protein (25-260AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arginine esterase (EC 3.4.21.35)
Immunogen Species	Dog
Research Area	Others

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: Arginine esterase antibody at 3μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 42 kDa

observed band size: 42 kDa