




zapA Antibody

Product Code	CSB-PA423441HA01BOL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8FG14
Immunogen	Recombinant Bacillus pumilus Cell division protein ZapA protein (1-85AA)
Raised In	Rabbit
Species Reactivity	Bacillus pumilus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:2000-1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division protein ZapA (Z ring-associated protein ZapA), zapA
Immunogen Species	Bacillus pumilus
Research Area	Others
Target Names	zapA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: zapA antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 26 kDa

Observed band size: 26 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.