



zapA Antibody, HRP conjugated

Product Code	CSB-PA423441HB01BOL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8FG14
Immunogen	Recombinant Bacillus pumilus Cell division protein ZapA protein (1-85AA)
Raised In	Rabbit
Species Reactivity	Bacillus pumilus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division protein ZapA (Z ring-associated protein ZapA), zapA
Immunogen Species	Bacillus pumilus (strain SAFR-032)
Research Area	Others
Target Names	zapA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.