




pyrE Antibody

Product Code	CSB-PA506933HA01LNO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C1D6F5
Immunogen	Recombinant Laribacter hongkongensis Orotate phosphoribosyltransferase protein (1-213AA)
Raised In	Rabbit
Species Reactivity	Laribacter hongkongensis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Orotate phosphoribosyltransferase (OPRT) (OPRTase) (EC 2.4.2.10), pyrE
Immunogen Species	Laribacter hongkongensis
Research Area	Others
Target Names	pyrE

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: pyrE antibody at 2.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 23 kDa

observed band size: 23 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.