




IplA Antibody

Product Code	CSB-PA510777HA01ENU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C4ZT68
Immunogen	Recombinant Escherichia coli Lipoate-protein ligase A protein (1-338AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipoate-protein ligase A (EC 6.3.1.20) (Lipoate--protein ligase), IplA
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	IplA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: IplA antibody at 2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 38 kDa

observed band size: 38 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.