



IplA Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA510777HB01ENU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | C4ZT68 |
| Immunogen | Recombinant Escherichia coli Lipoate-protein ligase A protein (1-338AA) |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Lipoate-protein ligase A (EC 6.3.1.20) (Lipoate--protein ligase), IplA |
| Immunogen Species | Escherichia coli (strain K12 / MC4100 / BW2952) |
| Research Area | Others |
| Target Names | IplA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |