



IplA Antibody, FITC conjugated

Product Code	CSB-PA510777HC01ENU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C4ZT68
Immunogen	Recombinant Escherichia coli Lipoate-protein ligase A protein (1-338AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipoate-protein ligase A (EC 6.3.1.20) (Lipoate--protein ligase), IplA
Immunogen Species	Escherichia coli (strain K12 / MC4100 / BW2952)
Research Area	Others
Target Names	IplA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.