




AK Antibody

Product Code	CSB-PA511311HA01ETF
Abbreviation	Arginine kinase
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C7E3T4
Immunogen	Recombinant Penaeus monodon Arginine kinase protein (2-356AA)
Raised In	Rabbit
Species Reactivity	Penaeus monodon
Tested Applications	ELISA, WB; Recommended dilution: WB:1:2000-1:10000
Relevance	arginine kinase activity, phosphorylation
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arginine kinase (AK) (EC 2.7.3.3) (allergen Pen m 2), AK
Species	Penaeus monodon
Research Area	Others
Target Names	AK

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: AK antibody at 2.8µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 55 kDa

Observed band size: 55 kDa