



AK Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA511311HC01ETF |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | C7E3T4 |
| Immunogen | Recombinant Penaeus monodon Arginine kinase protein (2-356AA) |
| Raised In | Rabbit |
| Species Reactivity | Penaeus monodon |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Arginine kinase (AK) (EC 2.7.3.3) (allergen Pen m 2), AK |
| Immunogen Species | Penaeus monodon (Giant tiger prawn) |
| Research Area | Others |
| Target Names | AK |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |