



tir Antibody

Product Code	CSB-PA514896HA01EOF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C6UYL8
Immunogen	Recombinant Escherichia coli O157:H7 Translocated intimin receptor Tir protein (252-362AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli O157:H7
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Translocated intimin receptor Tir (Secreted effector protein Tir), tir, espE
Immunogen Species	Escherichia coli O157:H7 (strain TW14359 / EHEC)
Research Area	Others
Target Names	tir
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.