



# PRO1 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA517332HB01EQB
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O04725
<b>Immunogen</b>	Recombinant Cynodon dactylon Profilin protein (2-131AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Cynodon dactylon
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Profilin (Pollen allergen Cyn d 12) (allergen Cyn d 12), PRO1; PRO2
<b>Immunogen Species</b>	Cynodon dactylon (Bermuda grass) (Panicum dactylon)
<b>Research Area</b>	Others
<b>Target Names</b>	PRO1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.