



PRO1 Antibody, FITC conjugated

Product Code	CSB-PA517332HC01EQB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O04725
Immunogen	Recombinant Cynodon dactylon Profilin protein (2-131AA)
Raised In	Rabbit
Species Reactivity	Cynodon dactylon
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Profilin (Pollen allergen Cyn d 12) (allergen Cyn d 12), PRO1; PRO2
Immunogen Species	Cynodon dactylon (Bermuda grass) (Panicum dactylon)
Research Area	Others
Target Names	PRO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.