




aspf3 Antibody

Product Code	CSB-PA525004HA01NGS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43099
Immunogen	Recombinant Neosartorya fumigata Peroxiredoxin Asp f3 protein (1-168AA)
Raised In	Rabbit
Species Reactivity	Neosartorya fumigata
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxiredoxin Asp f3 (Prx) (EC 1.11.1.15) (Thioredoxin peroxidase) (TPx) (allergen Asp f 3), aspf3
Immunogen Species	Neosartorya fumigata
Research Area	Others
Target Names	aspf3

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: aspf3 antibody at 3.25µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 35 kDa

Observed band size: 35 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.