



aspf3 Antibody, HRP conjugated

Product Code	CSB-PA525004HB01NGS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43099
Immunogen	Recombinant Neosartorya fumigata Peroxiredoxin Asp f3 protein (1-168AA)
Raised In	Rabbit
Species Reactivity	Neosartorya fumigata
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxiredoxin Asp f3 (Prx) (EC 1.11.1.15) (Thioredoxin peroxidase) (TPx) (allergen Asp f 3), aspf3
Immunogen Species	Neosartorya fumigata
Research Area	Others
Target Names	aspf3