



# aspf3 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA525004HC01NGS
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O43099
<b>Immunogen</b>	Recombinant Neosartorya fumigata Peroxiredoxin Asp f3 protein (1-168AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Neosartorya fumigata
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Peroxiredoxin Asp f3 (Prx) (EC 1.11.1.15) (Thioredoxin peroxidase) (TPx) (allergen Asp f 3), aspf3
<b>Immunogen Species</b>	Neosartorya fumigata
<b>Research Area</b>	Others
<b>Target Names</b>	aspf3
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.