



Profilin Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA525079HC01MQZ |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O49894 |
| Immunogen | Recombinant Mercurialis annua Profilin protein (1-133AA) |
| Raised In | Rabbit |
| Species Reactivity | Mercurialis annua |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Profilin (Pollen allergen Mer a 1) (allergen Mer a 1) |
| Immunogen Species | Mercurialis annua (Annual mercury) |
| Research Area | Others |