



NUPR1 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA527662LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O60356 |
| Immunogen | Recombinant Human Nuclear protein 1 protein (2-64AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Nuclear protein 1 (Candidate of metastasis 1) (Nuclear transcriptional regulator protein 1) (Protein p8), NUPR1, COM1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | NUPR1 |