



NUPR1 Antibody, FITC conjugated

Product Code	CSB-PA527662LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60356
Immunogen	Recombinant Human Nuclear protein 1 protein (2-64AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear protein 1 (Candidate of metastasis 1) (Nuclear transcriptional regulator protein 1) (Protein p8), NUPR1, COM1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NUPR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.