



NUPR1 Antibody, Biotin conjugated

Product Code	CSB-PA527662LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60356
Immunogen	Recombinant Human Nuclear protein 1 protein (2-64AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear protein 1 (Candidate of metastasis 1) (Nuclear transcriptional regulator protein 1) (Protein p8), NUPR1, COM1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NUPR1