



cobB Antibody, FITC conjugated

Product Code	CSB-PA536037LC01XAZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B0RM75
Immunogen	Recombinant Xanthomonas campestris pv. campestris NAD-dependent protein deacetylase protein (1-293AA)
Raised In	Rabbit
Species Reactivity	Xanthomonas campestris pv. campestris
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD-dependent protein deacetylase (EC 3.5.1.-) (Regulatory protein SIR2 homolog), cobB
Immunogen Species	Xanthomonas campestris pv. campestris
Research Area	Others
Target Names	cobB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.