




csrA Antibody

Product Code	CSB-PA543019HA01ENT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1XCM4
Immunogen	Recombinant Escherichia coli Translational regulator CsrA protein (1-61AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Translational regulator CsrA (Carbon storage regulator), csrA
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	csrA

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: csrA antibody at 3.4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 23 kDa

Observed band size: 23 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.