



GOLGA3 Antibody, HRP conjugated

Product Code	CSB-PA600096LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08378
Immunogen	Recombinant Human Golgin subfamily A member 3 protein (217-405AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Golgin subfamily A member 3 (Golgi complex-associated protein of 170 kDa) (GCP170) (Golgin-160), GOLGA3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GOLGA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.